

Wireless Parameters for Starkey Hearing Aids

Omega AI Family

Starkey wireless hearing aids contain a radio transceiver utilizing Bluetooth Low Energy wireless technology operating in the 2.4-2.4835 GHz frequency band with a maximum effective radiated power listed in the following table using GFSK transmission modulation. The receiver section of the radio has a bandwidth of 1.5 MHz.

Hearing aids with NFMI checked in the following table also contain a radio transceiver utilizing Near Field Magnetic Induction operating on 10.281 MHz with maximum induced magnetic field strength of -5 dBuA/m at a measurement distance of 10 meters with 8-DPSK transmission modulation. The receiver section of the NFMI radio has a bandwidth of 400 kHz.

Hearing aids with Telecoil checked are equipped with a telecoil. For custom ITE hearing aids, available means this hearing aid style can be equipped with a telecoil.

The EU Radio Equipment Directive Declaration of Conformity for each style and brand can be found at docs.starkeyhearingtechnologies.com.

A gray cell means that that the specified brand and style is not sold in the listed jurisdiction.

Style: RIC RT (with LED)										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Omega AI RIC RT	-1.5	√	√	EOA-24EDGRICRT	6903A-24EDGRICRT	Omega AI RIC RT	  205-240012	 R-R-EOA-AI24RICRT	 CCAB25LP0310T2
Audibel	Aris AI RIC RT	-1.5	√	√	EOA-24EDGRICRT	6903A-24EDGRICRT	Aris AI RIC RT	  205-240012	 R-R-EOA-AI24RICRT	 CCAB25LP031AT4
NuEar	Noble AI RIC RT	-1.5	√	√	EOA-24EDGRICRT	6903A-24EDGRICRT	Noble AI RIC RT	  205-240012	 R-R-EOA-AI24RICRT	 CCAB25LP031BT6

This applies to the version of the Omega AI RIC RT with the LED from the Chaska program.

Wireless Parameters for Starkey Hearing Aids

Omega AI Family

Style: mRIC R (with LED)										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Omega AI mRIC R	-5.6	X	X	EOA-24EDGMRIC	6903A-24EDGMRIC	Omega AI mRIC R	  205-240014	 R-R-EOA-AI24RICRT	 CCAB25LP0280T6
Audibel	Audibel Aris AI mRIC R	-5.6	X	X	EOA-24EDGMRIC	6903A-24EDGMRIC	Aris AI mRIC R	  205-240014	 R-R-EOA-AI24RICRT	 CCAB25LP028AT8
NuEar	NuEar Noble AI mRIC R	-5.6	X	X	EOA-24EDGMRIC	6903A-24EDGMRIC	Noble AI mRIC R	  205-240014	 R-R-EOA-AI24RICRT	 CCAB25LP028BT0

This applies to the version of the Omega AI mRIC R with the LED from the Chaska program.

Starkey, Starkey Edge, and StarLink are trademarks of Starkey Laboratories, Inc.

©2025 Starkey Laboratories, Inc. All Rights Reserved. 10/25 SPEC3094-02-EN-XX

Wireless Parameters for Starkey Hearing Aids

Omega AI Family

Style: RIC 312										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Omega AI RIC 312	-2.3	√	X	EOA-24EDGR312	6903A-24EDGR312	Omega AI RIC 312	  205-240013	 R-R-EOA-AI24RICRT	 CCAB25LP0320T5
Audibel	Audibel Aris AI RIC 312	-2.3	√	X	EOA-24EDGR312	6903A-24EDGR312	Aris AI RIC 312	  205-240013	 R-R-EOA-AI24RICRT	 CCAB25LP032AT7
NuEar	Nu Ear Noble AI RIC 312	-2.3	√	X	EOA-24EDGR312	6903A-24EDGR312	Noble AI RIC 312	  205-240013	 R-R-EOA-AI24RICRT	 CCAB25LP032BT9

Wireless Parameters for Starkey Hearing Aids

Omega AI Family

Style: ITE R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Omega AI ITE R	-1.7	X	Available	EOA-24EDGITELCH	6903A-24EDGITELCH	Omega AI ITE L	  205-240010	 R-R-SGk-24LVI0ITER	 CCAB25LP0290T9
Audibel	Audibel Aris AI ITE R	-1.7	X	Available	EOA-24EDGITELCH	6903A-24EDGITELCH	Aris AI ITE L	  205-240010	 R-R-SGk-24LVI0ITER	 CCAB25LP029AT1
NuEar	NuEar Noble AI ITE R	-1.7	X	Available	EOA-24EDGITELCH	6903A-24EDGITELCH	Noble AI ITE L	  205-240010	 R-R-SGk-24LVI0ITER	 CCAB25LP029BT0

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Starkey, Starkey Edge, and StarLink are trademarks of Starkey Laboratories, Inc.

©2025 Starkey Laboratories, Inc. All Rights Reserved. 10/25 SPEC3094-02-EN-XX

Wireless Parameters for Starkey Hearing Aids

Omega AI Family

Style: ITE R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Omega AI ITE R	-1.6	X	Available	EOA-24EDGITERCH	6903A-24EDGITERCH	Omega AI ITE R	  205-240011	 R-R-SGk-24LVI0ITER	 CCAB25LP0300T2
Audibel	Audibel Aris AI ITE R	-1.6	X	Available	EOA-24EDGITERCH	6903A-24EDGITERCH	Aris AI ITE R	  205-240011	 R-R-SGk-24LVI0ITER	 CCAB25LP030AT1
NuEar	NuEar Noble AI ITE R	-1.6	X	Available	EOA-24EDGITERCH	6903A-24EDGITERCH	Noble AI ITE R	  205-240011	 R-R-SGk-24LVI0ITER	 CCAB25LP030BT3

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Starkey, Starkey Edge, and StarLink are trademarks of Starkey Laboratories, Inc.

©2025 Starkey Laboratories, Inc. All Rights Reserved. 10/25 SPEC3094-02-EN-XX

Wireless Parameters for Starkey Hearing Aids

Omega AI Family

Style: ITC R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Omega AI ITC R	-1.7	X	X	EOA-24EDGITELCH	6903A-24EDGITELCH	Omega AI ITC L	  205-240010	 R-R-SGk-24LVM0ITER	 CCAB25LP0290T9
Audibel	Audibel Aris AI ITC R	-1.7	X	X	EOA-24EDGITELCH	6903A-24EDGITELCH	Aris AI ITC L	  205-240010	 R-R-SGk-24LVM0ITER	 CCAB25LP029AT1
NuEar	NuEar Noble AI ITC R	-1.7	X	X	EOA-24EDGITELCH	6903A-24EDGITELCH	Noble AI ITC L	  205-240010	 R-R-SGk-24LVM0ITER	 CCAB25LP029BT0

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Wireless Parameters for Starkey Hearing Aids

Omega AI Family

Style: ITC R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Omega AI ITC R	-1.6	X	X	EOA-24EDGITERCH	6903A-24EDGITERCH	Omega AI ITC R	  205-240011	 R-R-SGk-24LVI0ITER	 CCAB25LP0300T2
Audibel	Audibel Aris AI ITC R	-1.6	X	X	EOA-24EDGITERCH	6903A-24EDGITERCH	Aris AI ITC R	  205-240011	 R-R-SGk-24LVI0ITER	 CCAB25LP030AT1
NuEar	NuEar Noble AI ITC R	-1.6	X	X	EOA-24EDGITERCH	6903A-24EDGITERCH	Noble AI ITC R	  205-240011	 R-R-SGk-24LVI0ITER	 CCAB25LP030BT3

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Starkey, Starkey Edge, and StarLink are trademarks of Starkey Laboratories, Inc.

©2025 Starkey Laboratories, Inc. All Rights Reserved. 10/25 SPEC3094-02-EN-XX

Wireless Parameters for Starkey Hearing Aids

Omega AI Family

Style: CIC										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Omega AI CIC	0.4	X	X	EOA-24EDGCIC	6903A-24EDGCIC	Omega AI CIC	 205-240015	 R-R-SGk-24LVI0ITER	 CCAB25LP0330T8
Audibel	Audibel Aris AI CIC	0.4	X	X	EOA-24EDGCIC	6903A-24EDGCIC	Aris AI CIC	 205-240015	 R-R-SGk-24LVI0ITER	 CCAB25LP033AT0
NuEar	NuEar Noble AI CIC	0.4	X	X	EOA-24EDGCIC	6903A-24EDGCIC	Noble AI CIC	 205-240015	 R-R-SGk-24LVI0ITER	 CCAB25LP033BT2

Wireless Parameters for Starkey Hearing Aids

Edge AI Family

Starkey wireless hearing aids contain a radio transceiver utilizing Bluetooth Low Energy wireless technology operating in the 2.4-2.4835 GHz frequency band with a maximum effective radiated power listed in the following table using GFSK transmission modulation. The receiver section of the radio has a bandwidth of 1.5 MHz.

Hearing aids with NFMI checked in the following table also contain a radio transceiver utilizing Near Field Magnetic Induction operating on 10.281 MHz with maximum induced magnetic field strength of -5 dBuA/m at a measurement distance of 10 meters with 8-DPSK transmission modulation. The receiver section of the NFMI radio has a bandwidth of 400 kHz.

The EU Radio Equipment Directive Declarations of Conformity for these hearing aids can be found at docs.starkeyhearingtechnologies.com.

A gray cell means that that the specified brand and style is not sold in the listed jurisdiction.

Style: RIC RT										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Edge AI RIC RT	-1.5	√	√	EOA-24EDGRICRT	6903A-24EDGRICRT	Edge AI RIC RT	  205-240012	 R-R-EOA-AI24RICRT	 CCAB25LP0310T2
Audibel	Audibel Vitality AI RIC RT	-1.5	√	√	EOA-24EDGRICRT	6903A-24EDGRICRT	Vitality AI RIC RT	  205-240012	 R-R-EOA-AI24RICRT	 CCAB25LP031AT4
NuEar	NuEar Edge AI RIC RT	-1.5	√	√	EOA-24EDGRICRT	6903A-24EDGRICRT	Edge AI RIC RT	  205-240012	 R-R-EOA-AI24RICRT	 CCAB25LP031BT6
Audigy	AGXs Edge AI RIC RT	-1.5	√	√	EOA-24EDGRICRT	6903A-24EDGRICRT	Edge AI RIC RT	  205-240012	 R-R-EOA-AI24RICRT	
SoundGear	SoundGear AMP RIC RT	-1.5	√	√	EOA-24EDGRICRT	6903A-24EDGRICRT	SoundGear AMP RIC RT	  205-240012	 R-R-EOA-AI24RICRT	

Wireless Parameters for Starkey Hearing Aids

Edge AI Family

Style: mRIC R										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Edge AI mRIC R	-5.6	X	X	EOA-24EDGMRIC	6903A-24EDGMRIC	Edge AI RIC mRIC R	  205-240014	 R-R-EOA-AI24RICRT	 CCAB25LP0280T6
Audibel	Audibel Vitality AI mRIC R	-5.6	X	X	EOA-24EDGMRIC	6903A-24EDGMRIC	Vitality AI mRIC R	  205-240014	 R-R-EOA-AI24RICRT	 CCAB25LP028AT8
NuEar	NuEar Edge AI mRIC R	-5.6	X	X	EOA-24EDGMRIC	6903A-24EDGMRIC	Edge AI mRIC R	  205-240014	 R-R-EOA-AI24RICRT	 CCAB25LP028BT0
Audigy	AGXs Edge AI mRIC R	-5.6	X	X	EOA-24EDGMRIC	6903A-24EDGMRIC	Edge AI mRIC R	  205-240014	 R-R-EOA-AI24RICRT	
SoundGear	SoundGear AMP mRIC R	-5.6	X	X	EOA-24EDGMRIC	6903A-24EDGMRIC	SoundGear AMP mRIC R	  205-240014	 R-R-EOA-AI24RICRT	
KIND	KINDvaro	-5.6	X	X	EOA-24EDGMRIC			  205-240014	 R-R-EOA-AI24RICRT	

Wireless Parameters for Starkey Hearing Aids

Edge AI Family

Style: RIC 312										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Edge AI RIC 312	-2.3	√	X	EOA-24EDGR312	6903A-24EDGR312	Edge AI RIC 312	  205-240013	 R-R-EOA-AI24RICRT	 CCAB25LP0320T5
Audibel	Audibel Vitality AI RIC 312	-2.3	√	X	EOA-24EDGR312	6903A-24EDGR312	Vitality AI RIC 312	  205-240013	 R-R-EOA-AI24RICRT	 CCAB25LP032AT7
NuEar	Nu Ear Edge AI RIC 312	-2.3	√	X	EOA-24EDGR312	6903A-24EDGR312	Edge AI RIC 312	  205-240013	 R-R-EOA-AI24RICRT	 CCAB25LP032BT9
Audigy	AGXs Edge AI RIC 312	-2.3	√	X	EOA-24EDGR312	6903A-24EDGR312	Edge AI RIC 312	  205-240013	 R-R-EOA-AI24RICRT	
SoundGear	SoundGear AMP RIC 312	-2.3	√	X	EOA-24EDGR312	6903A-24EDGR312	SoundGear AMP RIC 312	  205-240013	 R-R-EOA-AI24RICRT	
KIND	KINDvaro	-2.3	√	X	EOA-24EDGR312			  205-240013	 R-R-EOA-AI24RICRT	

Wireless Parameters for Starkey Hearing Aids

Edge AI Family

Style: ITE R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Edge AI ITE R	-1.7	X	Available	EOA-24EDGITELCH	6903A-24EDGITELCH	Edge AI ITE L	  205-240010	 R-R-SGk-24LMOITER	 CCAB25LP0290T9
Audibel	Audibel Vitality AI ITE R	-1.7	X	Available	EOA-24EDGITELCH	6903A-24EDGITELCH	Vitality AI ITE L	  205-240010	 R-R-SGk-24LMOITER	 CCAB25LP029AT1
NuEar	NuEar Edge AI ITE R	-1.7	X	Available	EOA-24EDGITELCH	6903A-24EDGITELCH	Edge AI ITE L	  205-240010	 R-R-SGk-24LMOITER	 CCAB25LP029BT0
Audigy	AGXs Edge AI ITE R	-1.7	X	Available	EOA-24EDGITELCH	6903A-24EDGITELCH	Edge AI ITE L	  205-240010	 R-R-SGk-24LMOITER	
SoundGear	SoundGear AMP ITE R	-1.7	X	Available	EOA-24EDGITELCH	6903A-24EDGITELCH	SoundGear AMP ITE L	  205-240010	 R-R-SGk-24LMOITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Wireless Parameters for Starkey Hearing Aids

Edge AI Family

Style: ITE R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Edge AI ITE R	-1.6	X	Available	EOA-24EDGITERCH	6903A-24EDGITERCH	Edge AI ITE R	  205-240011	 R-R-SGk-24LMOITER	 CCAB25LP0300T2
Audibel	Audibel Vitality AI ITE R	-1.6	X	Available	EOA-24EDGITERCH	6903A-24EDGITERCH	Vitality AI ITE R	  205-240011	 R-R-SGk-24LMOITER	 CCAB25LP030AT1
NuEar	NuEar Edge AI ITE R	-1.6	X	Available	EOA-24EDGITERCH	6903A-24EDGITERCH	Edge AI ITE R	  205-240011	 R-R-SGk-24LMOITER	 CCAB25LP030BT3
Audigy	AGXs Edge AI ITE R	-1.6	X	Available	EOA-24EDGITERCH	6903A-24EDGITERCH	Edge AI ITE R	  205-240011	 R-R-SGk-24LMOITER	
SoundGear	SoundGear AMP ITE R	-1.6	X	Available	EOA-24EDGITERCH	6903A-24EDGITERCH	SoundGear AMP ITE R	  205-240011	 R-R-SGk-24LMOITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Wireless Parameters for Starkey Hearing Aids

Edge AI Family

Style: ITC R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Edge AI ITC R	-1.7	X	X	EOA-24EDGITECH	6903A-24EDGITECH	Edge AI ITC L	  205-240010	 R-R-SGk-24LMOITER	 CCAB25LP0290T9
Audibel	Audibel Vitality AI ITC R	-1.7	X	X	EOA-24EDGITECH	6903A-24EDGITECH	Vitality AI ITC L	  205-240010	 R-R-SGk-24LMOITER	 CCAB25LP029AT1
NuEar	NuEar Edge AI ITC R	-1.7	X	X	EOA-24EDGITECH	6903A-24EDGITECH	Edge AI ITC L	  205-240010	 R-R-SGk-24LMOITER	 CCAB25LP029BT0
Audigy	AGXs Edge AI ITC R	-1.7	X	X	EOA-24EDGITECH	6903A-24EDGITECH	Edge AI ITC L	  205-240010	 R-R-SGk-24LMOITER	
SoundGear	SoundGear AMP ITC R	-1.7	X	X	EOA-24EDGITECH	6903A-24EDGITECH	SoundGear AMP ITC L	  205-240010	 R-R-SGk-24LMOITER	
KIND	KINDvaro	-1.7	X	X	EOA-24EDGITECH			  205-240010	 R-R-SGk-24LMOITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Wireless Parameters for Starkey Hearing Aids

Edge AI Family

Style: ITC R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Edge AI ITC R	-1.6	X	X	EOA-24EDGITERCH	6903A-24EDGITERCH	Edge AI ITC R	  205-240011	 R-R-SGk-24LMOITER	 CCAB25LP0300T2
Audibel	Audibel Vitality AI ITC R	-1.6	X	X	EOA-24EDGITERCH	6903A-24EDGITERCH	Vitality AI ITC R	  205-240011	 R-R-SGk-24LMOITER	 CCAB25LP030AT1
NuEar	NuEar Edge AI ITC R	-1.6	X	X	EOA-24EDGITERCH	6903A-24EDGITERCH	Edge AI ITC R	  205-240011	 R-R-SGk-24LMOITER	 CCAB25LP030BT3
Audigy	AGXs Edge AI ITC R	-1.6	X	X	EOA-24EDGITERCH	6903A-24EDGITERCH	Edge AI ITC R	  205-240011	 R-R-SGk-24LMOITER	
SoundGear	SoundGear AMP ITC R	-1.6	X	X	EOA-24EDGITERCH	6903A-24EDGITERCH	SoundGear AMP ITC R	  205-240011	 R-R-SGk-24LMOITER	
KIND	KINDvaro	-1.6	X	X	EOA-24EDGITERCH			  205-240011	 R-R-SGk-24LMOITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Wireless Parameters for Starkey Hearing Aids

Edge AI Family

Style: CIC										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Edge AI CIC	0.4	X	X	EOA-24EDGCIC	6903A-24EDGCIC	Edge AI CIC	  205-240015	 R-R-SGk-24LMOITER	 CCAB25LP0330T8
Audibel	Audibel Vitality AI CIC	0.4	X	X	EOA-24EDGCIC	6903A-24EDGCIC	Vitality AI CIC	  205-240015	 R-R-SGk-24LMOITER	 CCAB25LP033AT0
NuEar	NuEar Edge AI CIC	0.4	X	X	EOA-24EDGCIC	6903A-24EDGCIC	Edge AI CIC	  205-240015	 R-R-SGk-24LMOITER	 CCAB25LP033BT2
Audigy	AGXs Edge AI CIC	0.4	X	X	EOA-24EDGCIC	6903A-24EDGCIC	Edge AI CIC	  205-240015	 R-R-SGk-24LMOITER	
SoundGear	SoundGear AMP CIC	0.4	X	X	EOA-24EDGCIC	6903A-24EDGCIC	SoundGear AMP CIC	  205-240015	 R-R-SGk-24LMOITER	
KIND	KINDvaro	0.4	X	X	EOA-24EDGCIC			  205-240015	 R-R-SGk-24LMOITER	

Wireless Parameters for Starkey Hearing Aids

Signature CIC R NW

Starkey wireless hearing aids contain a radio transceiver utilizing Bluetooth Low Energy wireless technology operating in the 2.4-2.4835 GHz frequency band with a maximum effective radiated power listed in the following table using GFSK transmission modulation. The receiver section of the radio has a bandwidth of 1.5 MHz.

The EU *Radio Equipment Directive* Declaration of Conformity for each style and brand can be found at docs.starkeyhearingtechnologies.com.

A gray cell means that that the specified brand and style is not sold in the listed jurisdiction.

Style: CIC R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Signature CIC R NW*	-8	X	X	EOA-24GENCICLPO	6903A-24GENCICLPO	Signature CIC R Left	  205-240001	 R-R-SGk-24LMOITER	 CCAB24LP0740T6
Audibel	Audibel Signature CIC R NW*	-8	X	X	EOA-24GENCICLPO	6903A-24GENCICLPO	Signature CIC R Left	  205-240001	 R-R-SGk-24LMOITER	 CCAB24LP074AT8
NuEar	NuEar Signature CIC R NW*	-8	X	X	EOA-24GENCICLPO	6903A-24GENCICLPO	Signature CIC R Left	  205-240001	 R-R-SGk-24LMOITER	 CCAB24LP074BT0
Audigy	AGXs Signature CIC R NW*	-8	X	X	EOA-24GENCICLPO	6903A-24GENCICLPO	Signature CIC R Left	  205-240001	 R-R-SGk-24LMOITER	
KIND	KINDsignature CIC R	-8	X	X	EOA-24GENCICLPO			  205-240001	 R-R-SGk-24LMOITER	

*There are separate regulatory IDs for the left and right hearing aids in the CIC R style.

**The CIC R NW is a wireless hearing aid for regulatory purposes since it contains a radio transmitter. The wireless functionality is only used during the programming of the hearing aid. Since the Signature CIC R NW is a wireless hearing aid, wireless regulatory IDs apply to it.

Wireless Parameters for Starkey Hearing Aids

Signature CIC R NW

Style: CIC R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Signature CIC R NW	-8	X	X	EOA-24GENCICRPO	6903A-24GENCICRPO	Signature CIC R Right	  205-240002	 R-R-SGk-24LMOITER	 CCAB24LP0750T9
Audibel	Audibel Signature CIC R NW	-8	X	X	EOA-24GENCICRPO	6903A-24GENCICRPO	Signature CIC R Right	  205-240002	 R-R-SGk-24LMOITER	 CCAB24LP075AT1
NuEar	NuEar Signature CIC R NW	-8	X	X	EOA-24GENCICRPO	6903A-24GENCICRPO	Signature CIC R Right	  205-240002	 R-R-SGk-24LMOITER	 CCAB24LP075BT0
Audigy	AGXs Signature CIC R NW	-8	X	X	EOA-24GENCICRPO	6903A-24GENCICRPO	Signature CIC R Right	  205-240002	 R-R-SGk-24LMOITER	
KIND	KINDsignature CIC R	-8	X	X	EOA-24GENCICRPO			  205-240002	 R-R-SGk-24LMOITER	

*There are separate regulatory IDs for the left and right hearing aids in the CIC R style.

**The CIC R NW is a wireless hearing aid for regulatory purposes since it contains a radio transmitter. The wireless functionality is only used during the programming of the hearing aid. Since the Signature CIC R NW is a wireless hearing aid, wireless regulatory IDs apply to it.

Starkey, Starkey Edge, and StarLink are trademarks of Starkey Laboratories, Inc.

©2025 Starkey Laboratories, Inc. All Rights Reserved. 10/25 SPEC3094-02-EN-XX

Wireless Parameters for Starkey Hearing Aids

Genesis AI Family

Starkey wireless hearing aids contain a radio transceiver utilizing Bluetooth Low Energy wireless technology operating in the 2.4-2.4835 GHz frequency band with a maximum effective radiated power listed in the following table using GFSK transmission modulation. The receiver section of the radio has a bandwidth of 1.5 MHz

Hearing aids with NFMI checked in the following table also contain a radio transceiver utilizing Near Field Magnetic Induction operating on 10.281 MHz with maximum induced magnetic field strength of -5 dBuA/m at a measurement distance of 10 meters with 8-DPSK transmission modulation. The receiver section of the NFMI radio has a bandwidth of 400 kHz.

The EU *Radio Equipment Directive* Declaration of Conformity for each style and brand can be found at docs.starkeyhearingtechnologies.com.

A gray cell means that that the specified brand and style is not sold in the listed jurisdiction.

Style: RIC RT										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Genesis AI RIC RT	-6.2	X	√	EOA-24GENMACRO	6903A-24GENMACRO	Genesis AI RIC RT	  205-220007	 R-R-EOA-AI24RICRT	 CCAB24LP0600T3
Audibel	Audibel Intrigue AI RIC RT	-6.2	X	√	EOA-24GENMACRO	6903A-24GENMACRO	Intrigue AI RIC RT	  205-220007	 R-R-EOA-AI24RICRT	 CCAB24LP060AT5
NuEar	NuEar NE NXG AI RIC RT	-6.2	X	√	EOA-24GENMACRO	6903A-24GENMACRO	NE NXG AI RIC RT	  205-220007	 R-R-EOA-AI24RICRT	 CCAB24LP060CT9
Audigy	AGXs Genesis AI RIC RT	-6.2	X	√	EOA-24GENMACRO	6903A-24GENMACRO	AGXs Genesis AI RIC RT	  205-220007	 R-R-EOA-AI24RICRT	
KIND	KINDvaro	-6.2	X	√	EOA-24GENMACRO			  205-220007	 R-R-EOA-AI24RICRT	
ProAkustik	iPro GenAI RIC RT	-6.2	X	√	EOA-24GENMACRO			  205-220007	 R-R-EOA-AI24RICRT	
Telefunken	vation RIC RT	-6.2	X	√	EOA-24GENMACRO			  205-220007	 R-R-EOA-AI24RICRT	

Wireless Parameters for Starkey Hearing Aids

Genesis AI Family

Style: mRIC R										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Genesis AI mRIC R	-12.7	X	X	EOA-24GENMICRO	6903A-24GENMICRO	Genesis AI mRIC R	  205-220006	 R-R-EOA-AI24RICRT	 CCAB24LP0590T0
Audibel	Audibel Intrigue AI mRIC R	-12.7	X	X	EOA-24GENMICRO	6903A-24GENMICRO	Intrigue AI mRIC R	  205-220006	 R-R-EOA-AI24RICRT	 CCAB24LP059AT2
NuEar	NuEar NE NXG AI mRIC R	-12.7	X	X	EOA-24GENMICRO	6903A-24GENMICRO	NE NXG AI mRIC R	  205-220006	 R-R-EOA-AI24RICRT	 CCAB24LP059CT6
Audigy	AGXs Genesis AI mRIC R	-12.7	X	X	EOA-24GENMICRO	6903A-24GENMICRO	AGXs Genesis AI mRIC R	  205-220006	 R-R-EOA-AI24RICRT	
KIND	KINDvaro	-12.7	X	X	EOA-24GENMICRO			  205-220006	 R-R-EOA-AI24RICRT	
ProAkustik	iPro GenAI mRIC R	-12.7	X	X	EOA-24GENMICRO			  205-220006	 R-R-EOA-AI24RICRT	
Telefunken	vation mRIC R	-12.7	X	X	EOA-24GENMICRO			  205-220006	 R-R-EOA-AI24RICRT	

Wireless Parameters for Starkey Hearing Aids

Genesis AI Family

Style: RIC 312										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Genesis AI RIC 312	-4.2	√	X	EOA-24GENR312	6903A-24GENR312	Genesis AI RIC 312	  205-230067	 R-R-EOA-AI24RICRT	 CCAB24LP0670T1
Audibel	Audibel Intrigue AI RIC 312	-4.2	√	X	EOA-24GENR312	6903A-24GENR312	Intrigue AI RIC 312	  205-230067	 R-R-EOA-AI24RICRT	 CCAB24LP067AT0
NuEar	Nu Ear NE NXG AI RIC 312	-4.2	√	X	EOA-24GENR312	6903A-24GENR312	NE NXG AI RIC 312	  205-230067	 R-R-EOA-AI24RICRT	 CCAB24LP067CT4
Audigy	AGXs Genesis AI RIC 312	-4.2	√	X	EOA-24GENR312	6903A-24GENR312	AGXs Genesis AI RIC 312	  205-230067	 R-R-EOA-AI24RICRT	
KIND	KINDvaro	-4.2	√	X	EOA-24GENR312			  205-230067	 R-R-EOA-AI24RICRT	
ProAkustik	iPro GenAI RIC 312	-4.2	√	X	EOA-24GENR312			  205-230067	 R-R-EOA-AI24RICRT	
Telefunken	vation RIC 312	-4.2	√	X	EOA-24GENR312			  205-230067	 R-R-EOA-AI24RICRT	

Wireless Parameters for Starkey Hearing Aids

Genesis AI Family

Style: ITE R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Genesis AI ITE R	0	X	Available	EOA-24GENITELCH	6903A-24GENITELCH	Genesis AI ITE L	  205-240065	 R-R-SGk-24LVI0ITER	 CCAB24LP0720T0
Audibel	Audibel Intrigue AI ITE R	0	X	Available	EOA-24GENITELCH	6903A-24GENITELCH	Intrigue AI ITE L	  205-240065	 R-R-SGk-24LVI0ITER	 CCAB24LP072AT2
NuEar	NuEar NE NXG AI ITE R	0	X	Available	EOA-24GENITELCH	6903A-24GENITELCH	NE NXG AI ITE L	  205-240065	 R-R-SGk-24LVI0ITER	 CCAB24LP072CT6
Audigy	AGXs Genesis AI ITE R	0	X	Available	EOA-24GENITELCH	6903A-24GENITELCH	AGXs Genesis AI ITE L	  205-240065	 R-R-SGk-24LVI0ITER	
ProAkustik	iPro GenAI ITE R	0	X	Available	EOA-24GENITELCH			  205-240065	 R-R-SGk-24LVI0ITER	
Telefunken	vation ITE R	0	X	Available	EOA-24GENITELCH			  205-240065	 R-R-SGk-24LVI0ITER	

Style: ITE R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Genesis AI ITE R	0	X	Available	6903A-24GENITERCH	6903A-24GENITERCH	Genesis AI ITE R	  205-240066	 R-R-SGk-24LVI0ITER	 CCAB24LP0730T3
Audibel	Audibel Intrigue AI ITE R	0	X	Available	6903A-24GENITERCH	6903A-24GENITERCH	Intrigue AI ITE R	  205-240066	 R-R-SGk-24LVI0ITER	 CCAB24LP073AT5
NuEar	NuEar NE NXG AI ITE R	0	X	Available	6903A-24GENITERCH	6903A-24GENITERCH	NE NXG AI ITE R	  205-240066	 R-R-SGk-24LVI0ITER	 CCAB24LP073CT9
Audigy	AGXs Genesis AI ITE R	0	X	Available	6903A-24GENITERCH	6903A-24GENITERCH	AGXs Genesis AI ITE R	  205-240066	 R-R-SGk-24LVI0ITER	
ProAkustik	iPro GenAI ITE R	0	X	Available	6903A-24GENITERCH			  205-240066	 R-R-SGk-24LVI0ITER	
Telefunken	vation ITE R	0	X	Available	6903A-24GENITERCH			  205-240066	 R-R-SGk-24LVI0ITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles

Wireless Parameters for Starkey Hearing Aids

Genesis AI Family

Style: ITC R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Genesis AI ITC R	0	X	X	EOA-24GENITELCH	6903A-24GENITELCH	Genesis AI ITC L	  205-230065	 R-R-SGk-24LVI0ITER	 CCAB24LP0720T0
Audibel	Audibel Intrigue AI ITC R	0	X	X	EOA-24GENITELCH	6903A-24GENITELCH	Intrigue AI ITC L	  205-230065	 R-R-SGk-24LVI0ITER	 CCAB24LP072AT2
NuEar	NuEar NE NXG AI ITC R	0	X	X	EOA-24GENITELCH	6903A-24GENITELCH	NE NXG AI ITC L	  205-230065	 R-R-SGk-24LVI0ITER	 CCAB24LP072CT6
Audigy	AGXs Genesis AI ITC R	0	X	X	EOA-24GENITELCH	6903A-24GENITELCH	AGX Genesis AI ITC L	  205-230065	 R-R-SGk-24LVI0ITER	
KIND	KINDvaro	0	X	X	EOA-24GENITELCH			  205-230065	 R-R-SGk-24LVI0ITER	
ProAkustik	iPro GenAI ITC R	0	X	X	EOA-24GENITELCH			  205-230065	 R-R-SGk-24LVI0ITER	
Telefunken	vation ITC R	0	X	X	EOA-24GENITELCH			  205-230065	 R-R-SGk-24LVI0ITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Wireless Parameters for Starkey Hearing Aids

Genesis AI Family

Style: ITC R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Genesis AI ITC R	0	X	X	EOA-24GENITERCH	6903A-24GENITERCH	Genesis AI ITC R	  205-240066	 R-R-SGk-24LVI0ITER	 CCAB24LP0730T3
Audibel	Audibel Intrigue AI ITC R	0	X	X	EOA-24GENITERCH	6903A-24GENITERCH	Intrigue AI ITC R	  205-240066	 R-R-SGk-24LVI0ITER	 CCAB24LP073AT5
NuEar	NuEar NE NXG AI ITC R	0	X	X	EOA-24GENITERCH	6903A-24GENITERCH	NE NXG AI ITC R	  205-240066	 R-R-SGk-24LVI0ITER	 CCAB24LP073CT9
Audigy	AGXs Genesis AI ITC R	0	X	X	EOA-24GENITERCH	6903A-24GENITERCH	AGXs Genesis AI ITC R	  205-240066	 R-R-SGk-24LVI0ITER	
KIND	KINDvaro	0	X	X	EOA-24GENITERCH			  205-240066	 R-R-SGk-24LVI0ITER	
ProAkustik	iPro GenAI ITC R	0	X	X	EOA-24GENITERCH			  205-240066	 R-R-SGk-24LVI0ITER	
Telefunken	vation ITC R	0	X	X	EOA-24GENITERCH			  205-240066	 R-R-SGk-24LVI0ITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles.

Wireless Parameters for Starkey Hearing Aids

Genesis AI Family

Style: CIC										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Genesis AI CIC	0	X	X	EOA-24GENCICZN	6903A-24GENCICZN	Genesis AI CIC	 R 205-230064	 R-R-SGk-24LVI0ITER	 CCAB24LP0710T0
Audibel	Audibel Intrigue AI CIC	0	X	X	EOA-24GENCICZN	6903A-24GENCICZN	Intrigue AI CIC	 R 205-230064	 R-R-SGk-24LVI0ITER	 CCAB24LP071AT2
NuEar	NuEar NE NXG AI CIC	0	X	X	EOA-24GENCICZN	6903A-24GENCICZN	NE NXG AI CIC	 R 205-230064	 R-R-SGk-24LVI0ITER	 CCAB24LP071CT3
Audigy	AGXs Genesis AI CIC	0	X	X	EOA-24GENCICZN	6903A-24GENCICZN	AGXs Genesis AI CIC	 R 205-230064	 R-R-SGk-24LVI0ITER	
KIND	KINDvaro	0	X	X	EOA-24GENCICZN			 R 205-230064	 R-R-SGk-24LVI0ITER	
ProAkustik	iPro GenAI CIC	0	X	X	EOA-24GENCICZN			 R 205-230064	 R-R-SGk-24LVI0ITER	
Telefunken	vation CIC	0	X	X	EOA-24GENCICZN			 R 205-230064	 R-R-SGk-24LVI0ITER	

Wireless Parameters for Starkey Hearing Aids

G Series AI Family

Starkey wireless hearing aids contain a radio transceiver utilizing Bluetooth Low Energy wireless technology operating in the 2.4-2.4835 GHz frequency band with a maximum effective radiated power listed in the following table using GFSK transmission modulation. The receiver section of the radio has a bandwidth of 1.5 MHz.

Hearing aids with Telecoil checked are equipped with a telecoil.

The EU *Radio Equipment Directive* Declaration of Conformity for each style and brand can be found at docs.starkeyhearingtechnologies.com.

A gray cell means that that the specified brand and style is not sold in the listed jurisdiction.

Style: RIC R										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	G Series AI RIC R	-6.7	X	X	EOA-24GSERRICR	6903A-24GSERRICR	G Series AI RIC R	  205-250005	 R-R-EOA-AI24RICRT	 CCAB24LP0600T3
Audibel	A Series AI RIC R	-6.7	X	X	EOA-24GSERRICR	6903A-24GSERRICR	A Series AI RIC R	  205-250005	 R-R-EOA-AI24RICRT	 CCAB24LP060AT5
NuEar	NE Series AI RIC R	-6.7	X	X	EOA-24GSERRICR	6903A-24GSERRICR	NE Series AI RIC R	  205-250005	 R-R-EOA-AI24RICRT	 CCAB24LP060CT9

Wireless Parameters for Starkey Hearing Aids

G Series AI Family

Style: RIC 312										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	G Series AI RIC 312	-1.3	X	X	EOA-24GSERR312	6903A-24GSERR312	G Series AI RIC 312	 205-250006	 R-R-EOA-AI24RICRT	 CCAB24LP0670T1
Audibel	A Series AI RIC 312	-1.3	X	X	EOA-24GSERR312	6903A-24GSERR312	A Series AI RIC 312	 205-250006	 R-R-EOA-AI24RICRT	 CCAB24LP067AT0
NuEar	NE Series AI RIC 312	-1.3	X	X	EOA-24GSERR312	6903A-24GSERR312	NE Series AI RIC 312	 205-250006	 R-R-EOA-AI24RICRT	 CCAB24LP067CT4

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Starkey wireless hearing aids contain a radio transceiver utilizing Bluetooth Low Energy wireless technology operating in the 2.4-2.4835 GHz frequency band with a maximum effective radiated power listed in the following table using GFSK transmission modulation. The receiver section of the radio has a bandwidth of 1.5 MHz

Hearing aids with NFMI checked in the following table also contain a radio transceiver utilizing Near Field Magnetic Induction operating on 10.281 MHz with maximum induced magnetic field strength of -5 dBuA/m at a measurement distance of 10 meters with 8-DPSK transmission modulation. The receiver section of the NFMI radio has a bandwidth of 400 kHz.

The Grand Avenue document number for the EU *Radio Equipment Directive* Declaration of Conformity is shown for each style and brand. These Declarations of Conformity can also be found at docs.starkeyhearingtechnologies.com.

A gray cell means that that the specified brand and style is not sold in the listed jurisdiction.

Style: RIC 312										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI RIC 312	-13	√	X	EOA-24LMOR312	6903A-24LMOR312	Evolv AI RIC 312	  209-J00384 		
	Starkey Model S	-13	√	X	EOA-24LMOR312	6903A-24LMOR312	Model S RIC 312	  209-J00384 		
Audibel	Audibel Arc AI RIC 312	-13	√	X	EOA-24LMOR312	6903A-24LMOR312	Arc AI RIC 312	  209-J00384 		
	Audibel Model A	-13	√	X	EOA-24LMOR312	6903A-24LMOR312	Model A RIC 312	  209-J00384 		
NuEar	NuEar Savant AI RIC 312	-13	√	X	EOA-24LMOR312	6903A-24LMOR312	Savant AI RIC 312	  209-J00384 		
	NuEar Model N	-13	√	X	EOA-24LMOR312	6903A-24LMOR312	Model N RIC 312	  209-J00384 		
Audigy	AGXs Evolv AI RIC 312	-13	√	X	EOA-24LMOR312	6903A-24LMOR312	AGXs Evolv AI RIC 312	  209-J00384 		
	AGXss RIC 312	-13	√	X	EOA-24LMOR312	6903A-24LMOR312	Model AGXss RIC 312	  209-J00384 		
KIND	KINDvitalo	-13	√	X	EOA-24LMOR312			  209-J00384 		
ProAkustik	iPro 11.7 AI RIC 312	-13	√	X	EOA-24LMOR312			  209-J00384 		
	Model S iPro 1.5 RIC 312	-13	√	X	EOA-24LMOR312			  209-J00384 		
Telefunken	TF10-2+X RIC 312	-13	√	X	EOA-24LMOR312			  209-J00384 		
	Model S TF1-2 RIC 312	-13	√	X	EOA-24LMOR312			  209-J00384 		

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: mRIC 312										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI mRIC 312	-13	X	X	EOA-24LIVOM312	6903A-24LIVOM312	Evolv AI mRIC 312	  209-J00385	 MSIP-CRIM-EOA-24HAL02R312	 CCAM19LP1310T0
Audibel	Audibel Arc AI mRIC 312	-13	X	X	EOA-24LIVOM312	6903A-24LIVOM312	Arc AI RIC mRIC 312	  209-J00385	 MSIP-CRIM-EOA-24HAL02R312	 CCAM19LP1311T2
NuEar	NuEar Savant AI mRIC 312	-13	X	X	EOA-24LIVOM312	6903A-24LIVOM312	Savant AI mRIC 312	  209-J00385	 MSIP-CRIM-EOA-24HAL02R312	 CCAM19LP1313T3
Audigy	AGXs Evolv AI mRIC 312	-13	X	X	EOA-24LIVOM312	6903A-24LIVOM312	AGXs Evolv AI mRIC 312	  209-J00385	 MSIP-CRIM-EOA-24HAL02R312	
ProAkustik	iPro 11.7 AI mRIC 312	-13	X	X	EOA-24LIVOM312			  209-J00385	 MSIP-CRIM-EOA-24HAL02R312	
Telefunken	TF10-2+X mRIC 312	-13	X	X	EOA-24LIVOM312			  209-J00385	 MSIP-CRIM-EOA-24HAL02R312	

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: RIC R										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI RIC R	-13	√	√	EOA-24LMORCHG	6903A-24LMORCHG	Evolv AI RIC R	  209-J00383	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1321T2
Audibel	Audibel Arc AI RIC R	-13	√	√	EOA-24LMORCHG	6903A-24LMORCHG	Arc AI RIC R	  209-J00383	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1325T0
NuEar	Nu Ear Savant AI RIC R	-13	√	√	EOA-24LMORCHG	6903A-24LMORCHG	Savant AI RIC R	  209-J00383	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1327T1
Audigy	AGXs Evolv AI RIC R	-13	√	√	EOA-24LMORCHG	6903A-24LMORCHG	AGXS Evolv AI RIC R	  209-J00383	 MSIP-CRIM-EOA-24HAL02R13	
KIND	KINDvitalo	-13	√	√	EOA-24LMORCHG			  209-J00383	 MSIP-CRIM-EOA-24HAL02R13	
ProAkustik	iPro 11.7 AI RIC R	-13	√	√	EOA-24LMORCHG			  209-J00383	 MSIP-CRIM-EOA-24HAL02R13	
Telefunken	TF10-2+X RIC R	-13	√	√	EOA-24LMORCHG			  209-J00383	 MSIP-CRIM-EOA-24HAL02R13	

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: BTE 13										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI BTE 13	-9	✓	✓	EOA-24LMOB13A	6903A-24LMOB13A	Evolv AI *BTE13	 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1300T7
	Starkey Model S	-9	✓	✓	EOA-24LMOB13A	6903A-24LMOB13A	Model S BTE 13	 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
Audibel	Audibel Arc AI BTE 13	-9	✓	✓	EOA-24LMOB13A	6903A-24LMOB13A	Arc AI *BTE13	 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1301T9
	Audibel Model A	-9	✓	✓	EOA-24LMOB13A	6903A-24LMOB13A	Model A BTE 13	 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
NuEar	NuEar Savant AI BTE 13	-9	✓	✓	EOA-24LMOB13A	6903A-24LMOB13A	Savant AI *BTE13	 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1303T0
	NuEar Model N	-9	✓	✓	EOA-24LMOB13A	6903A-24LMOB13A	Model N BTE 13	 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
Audigy	AGXs Evolv AI BTE 13	-9	✓	✓	EOA-24LMOB13A	6903A-24LMOB13A	AGXs Evolv AI *BTE13	 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
	AGXss BTE 13	-9	✓	✓	EOA-24LMOB13A	6903A-24LMOB13A	Model AGXss *BTE13	 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
KIND	KINDvitalo	-9	✓	✓	EOA-24LMOB13A			 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
ProAkustik	iPro 11.7 AI RIC 312	-9	✓	✓	EOA-24LMOB13A			 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
	Model S iPro 1.5 RIC 312	-9	✓	✓	EOA-24LMOB13A			 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
Telefunken	TF10-2+X RIC 312	-9	✓	✓	EOA-24LMOB13A			 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	
	Model S TF1-2 RIC 312	-9	✓	✓	EOA-24LMOB13A			 R 209-J00387	 MSIP-CRIM-EOA-24HAL02R13	

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: BTE R										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI BTE R	-13	X	X	EOA-24LMI0BCHG	6903A-24LMI0BCHG	Evolv AI BTE R	 R 209-J00386	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1700T0
Audibel	Audibel Arc AI BTE R	-13	X	√	EOA-24LMI0BCHG	6903A-24LMI0BCHG	Arc AI BTE R	 R 209-J00386	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1701T2
NuEar	Nu Ear Savant AI BTE R	-13	X	√	EOA-24LMI0BCHG	6903A-24LMI0BCHG	Savant AI BTE R	 R 209-J00386	 MSIP-CRIM-EOA-24HAL02R13	 CCAM19LP1703T3
Audigy	AGXs Evolv AI BTE R	-13	X	√	EOA-24LMI0BCHG	6903A-24LMI0BCHG	AGXS Evolv AI BTE R	 R 209-J00386	 MSIP-CRIM-EOA-24HAL02R13	
KIND	KINDvitalo	-13	X	√	EOA-24LMI0BCHG			 R 209-J00386	 MSIP-CRIM-EOA-24HAL02R13	
ProAkustik	iPro 11.7 AI BTE R	-13	X	√	EOA-24LMI0BCHG			 R 209-J00386	 MSIP-CRIM-EOA-24HAL02R13	
Telefunken	TF10-2+X BTE R	-13	X	√	EOA-24LMI0BCHG			 R 209-J00386	 MSIP-CRIM-EOA-24HAL02R13	

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: Power Plus BTE 13										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI Power Plus BTE 13	-13	X	X	EOA-24LMOPB13	6903A-24LMOPB13	Evolv AI P+ BTE 13	  205-210103	 MSIP-CRIM-EOA-24HAL02R13	 CCAL21LP1170T7
Audibel	Audibel Arc AI Power Plus BTE 13	-13	X	√	EOA-24LMOPB13	6903A-24LMOPB13	Arc AI P+ BTE 13	  205-210103	 MSIP-CRIM-EOA-24HAL02R13	 CCAL21LP1171T9
NuEar	Nu Ear Savant AI Power Plus BTE 13	-13	X	√	EOA-24LMOPB13	6903A-24LMOPB13	Savant AI P+ BTE 13	  205-210103	 MSIP-CRIM-EOA-24HAL02R13	 CCAL21LP1173T0
Audigy	AGXs Evolv AI Power Plus BTE 13	-13	X	√	EOA-24LMOPB13	6903A-24LMOPB13	AGXS Evolv AI P+ BTE 13	  205-210103	 MSIP-CRIM-EOA-24HAL02R13	
KIND	KINDvitalo	-13	X	√	EOA-24LMOPB13			  205-210103	 MSIP-CRIM-EOA-24HAL02R13	
ProAkustik	iPro 11.7 AI P+ BTE 13	-13	X	√	EOA-24LMOPB13			  205-210103	 MSIP-CRIM-EOA-24HAL02R13	
Telefunken	TF10-2+X P+ BTE 13	-13	X	√	EOA-24LMOPB13			  205-210103	 MSIP-CRIM-EOA-24HAL02R13	

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: ITE R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI ITE R	-4	X	Available	EOA-24LVI0ITEL	6903A-24LVI0ITEL	Evolv AI ITEL	 R 209-J00403	 R-R-SGk-24LVI0ITER	 CCAM20LP0510T2
Audibel	Audibel Arc AI ITE R	-4	X	Available	EOA-24LVI0ITEL	6903A-24LVI0ITEL	Arc AI ITEL	 R 209-J00403	 R-R-SGk-24LVI0ITER	 CCAM20LP0511T4
NuEar	Nu Ear Savant AI ITE R	-4	X	Available	EOA-24LVI0ITEL	6903A-24LVI0ITEL	Savant AI ITEL	 R 209-J00403	 R-R-SGk-24LVI0ITER	 CCAM20LP0513T8
Audigy	AGXs Evolv AI ITE R	-4	X	Available	EOA-24LVI0ITEL	6903A-24LVI0ITEL	AGXS Evolv AI ITEL	 R 209-J00403	 R-R-SGk-24LVI0ITER	
ProAkustik	iPro 11.7 AI ITE R	-4	X	Available	EOA-24LVI0ITEL			 R 209-J00403	 R-R-SGk-24LVI0ITER	
Telefunken	TF10-2+X ITE R	-4	X	Available	EOA-24LVI0ITEL			 R 209-J00403	 R-R-SGk-24LVI0ITER	

Style: ITE R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI ITE R	-4	X	Available	EOA-24LVI0ITER	6903A-24LVI0ITER	Evolv AI ITER	 R 209-J00402	 R-R-SGk-24LVI0ITER	 CCAM20LP0520T5
Audibel	Audibel Arc AI ITE R	-4	X	Available	EOA-24LVI0ITER	6903A-24LVI0ITER	Arc AI ITER	 R 209-J00402	 R-R-SGk-24LVI0ITER	 CCAM20LP0521T7
NuEar	Nu Ear Savant AI ITE R	-4	X	Available	EOA-24LVI0ITER	6903A-24LVI0ITER	Savant AI ITER	 R 209-J00402	 R-R-SGk-24LVI0ITER	 CCAM20LP0523T1
Audigy	AGXs Evolv AI ITE R	-4	X	Available	EOA-24LVI0ITER	6903A-24LVI0ITER	AGXS Evolv AI ITER	 R 209-J00402	 R-R-SGk-24LVI0ITER	
ProAkustik	iPro 11.7 AI ITE R	-4	X	Available	EOA-24LVI0ITER			 R 209-J00402	 R-R-SGk-24LVI0ITER	
Telefunken	TF10-2+X ITE R	-4	X	Available	EOA-24LVI0ITER			 R 209-J00402	 R-R-SGk-24LVI0ITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: ITC R left*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI ITC R	-4	X	X	EOA-24LVIOITEL	6903A-24LVIOITEL	Evolv AI ITCL	  209-J00403	 R-R-SGk-24LVIOITER	 CCAM20LP0510T2
Audibel	Audibel Arc AI ITC R	-4	X	X	EOA-24LVIOITEL	6903A-24LVIOITEL	Arc AI ITCL	  209-J00403	 R-R-SGk-24LVIOITER	 CCAM20LP0511T4
NuEar	Nu Ear Savant AI ITC R	-4	X	X	EOA-24LVIOITEL	6903A-24LVIOITEL	Savant AI ITCL	  209-J00403	 R-R-SGk-24LVIOITER	 CCAM20LP0513T8
Audigy	AGXs Evolv AI ITC R	-4	X	X	EOA-24LVIOITEL	6903A-24LVIOITEL	AGXS Evolv AI ITCL	  209-J00403	 R-R-SGk-24LVIOITER	
KIND	KINDvitalo	-4	X	X	EOA-24LVIOITEL			  209-J00403	 R-R-SGk-24LVIOITER	
ProAkustik	iPro 11.7 AI ITC R	-4	X	X	EOA-24LVIOITEL			  209-J00403	 R-R-SGk-24LVIOITER	
Telefunken	TF10-2+X ITC R	-4	X	X	EOA-24LVIOITEL			  209-J00403	 R-R-SGk-24LVIOITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles

Starkey, Starkey Edge, and StarLink are trademarks of Starkey Laboratories, Inc.

©2025 Starkey Laboratories, Inc. All Rights Reserved. 10/25 SPEC3094-02-EN-XX

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: ITC R right*										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI ITC R	-4	X	X	EOA-24LIVIOITER	6903A-24LIVIOITER	Evolv AI ITCR	  209-J00402	 R-R-SGk-24LIVIOITER	 CCAM20LP0520T5
Audibel	Audibel Arc AI ITC R	-4	X	X	EOA-24LIVIOITER	6903A-24LIVIOITER	Arc AI ITCR	  209-J00402	 R-R-SGk-24LIVIOITER	 CCAM20LP0521T7
NuEar	Nu Ear Savant AI ITC R	-4	X	X	EOA-24LIVIOITER	6903A-24LIVIOITER	Savant AI ITCR	  209-J00402	 R-R-SGk-24LIVIOITER	 CCAM20LP0523T1
Audigy	AGXs Evolv AI ITC R	-4	X	X	EOA-24LIVIOITER	6903A-24LIVIOITER	AGXS Evolv AI ITCR	  209-J00402	 R-R-SGk-24LIVIOITER	
KIND	KINDvitalo	-4	X	X	EOA-24LIVIOITER			  209-J00402	 R-R-SGk-24LIVIOITER	
ProAkustik	iPro 11.7 AI ITC R	-4	X	X	EOA-24LIVIOITER			  209-J00402	 R-R-SGk-24LIVIOITER	
Telefunken	TF10-2+X ITC R	-4	X	X	EOA-24LIVIOITER			  209-J00402	 R-R-SGk-24LIVIOITER	

*There are separate regulatory IDs for the left and right hearing aids in the ITE and ITC styles

Starkey, Starkey Edge, and StarLink are trademarks of Starkey Laboratories, Inc.

©2025 Starkey Laboratories, Inc. All Rights Reserved. 10/25 SPEC3094-02-EN-XX

Wireless Parameters for Starkey Hearing Aids

Evolv AI Family

Style: CIC 312										
Brand	Brand Name	BT RF Power (dBm e.i.r.p)	NFMI	Telecoil	FCC ID	ISED IC	ISED HVIN	Japan ID	Korea ID	Taiwan ID
Starkey	Evolv AI CIC	0	X	X	EOA-24LIVIOCIC	6903A-24LIVIOCIC	Evolv AI CIC	  205-210102	 R-R-SGk-24LIVIOITER	 CCAL21LP1160T4
Audibel	Audibel Arc AI CIC	0	X	X	EOA-24LIVIOCIC	6903A-24LIVIOCIC	Arc AI CIC	  205-210102	 R-R-SGk-24LIVIOITER	 CCAL21LP1161T6
NuEar	Nu Ear Savant AI CIC	0	X	X	EOA-24LIVIOCIC	6903A-24LIVIOCIC	Savant AI CIC	  205-210102	 R-R-SGk-24LIVIOITER	 CCAL21LP1163T0
Audigy	AGXs Evolv AI CIC	0	X	X	EOA-24LIVIOCIC	6903A-24LIVIOCIC	AGXS Evolv AI CIC	  205-210102	 R-R-SGk-24LIVIOITER	
KIND	KINDvitalo	0	X	X	EOA-24LIVIOCIC			  205-210102	 R-R-SGk-24LIVIOITER	
ProAkustik	iPro 11.7 AI CIC	0	X	X	EOA-24LIVIOCIC			  205-210102	 R-R-SGk-24LIVIOITER	
Telefunken	TF10-2+X CIC	0	X	X	EOA-24LIVIOCIC			  205-210102	 R-R-SGk-24LIVIOITER	



Manufacturer according to FDA and EU Medical Device Regulation 2017/745

Starkey Laboratories, Inc.
6700 Washington Ave. South
Eden Prairie, MN 55344 USA
www.starkey.com

EC

REP

Starkey Laboratories (Germany) GmbH
Weg beim Jäger 218-222
22335 Hamburg
Germany

UK Authorised Person:
Starkey UK
William F. Austin House
Pepper Rd, Hazel Grove
Stockport SK7 5BX, UK
www.starkey.co.uk

Manufacturer according to Health Canada
Starkey Labs Canada Co.
2476 Argentia Road, Suite 301
Mississauga, ON L5N 6M1
www.starkeycanada.ca

CH

REP

QUNIQUE GmbH
Bahnhofweg 17
5610 Wohlen
Switzerland


2797

**Only applies to products sold in the EU*